Product Overview

Read the manual carefully before use, and return for future reference.

1. Mouth Protection Cover
2. On/Off Indicator
3. Brush Head Cover
4. Mouth Indicator
5. Charging Lock
6. Custom Mode

Do not share the brush head with others.
 Do not use the power cord as a means of charging dock while it is plugged in.
 Do not place any metal objects on top of the charging dock.
 Do not use the power cord as a means of charging dock while it is plugged in.
 Do not use the power cord as a means of charging dock while it is plugged in.
 Do not use the power cord as a means of charging dock while it is plugged in.

How to Use

1. Turn on / Turning On / Starting / Switching modes

2. Healthy brushing

3. Smart Functions

4. Smart Functions

5. Smart Functions

6. Smart Functions

Smart Functions

1. The toothbrush is equipped with an ultrasonic protection system that detects power in real time. When using the toothbrush, the power cord is automatically disconnected and the indicator will quickly flash red and the toothbrush will automatically power off.

2. You can set the brushing duration, modes and splash protection for the MI Home / Xiaomi Home Hub. After brushing with the MI Home / Xiaomi Home app, you can view the brushing data and generate a report.

3. Connecting with MI Home / Xiaomi Home app

4. The product is compatible with MI Home / Xiaomi Home app.

Specifications

Name: MI Smart Series Toothbrush T500 Model: M56501
 Brand Voltage: 3.7V ~ 4.2V Rated Power: 2W
 Rated Input: 5V 1.5A Charging Time: 12 hours
 Battery Capacity: 1000mAh
 Frequency Band: 2400-2480 MHz IP Rating: IPX7
 Battery Type: Rechargeable Lithium-Ion
 Dimensions: 26.1 x 24.6 x 85.0 mm

Charging

Item	Value
Charging Cable	300% (20%)
Charging Time	60% (20%)
Charging Speed	100%
Charging Lock	100% (20%)
Charging Lock	100% (20%)
Charging Lock	100% (20%)

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